


# AROLEAP PRO

The Aroleap Pro is a compact, high-precision device for rehabilitation and functional training. It provides adjustable assistance or resistance to support mobility, strength, and neuromuscular recovery. Therapists can apply controlled loading for safe, progressive rehab and therapeutic exercises. Real-time data tracks movement quality, range of motion, and load for accurate progress monitoring. Its versatile design fits easily in clinics, gyms, or home-care environments. The system helps improve recovery outcomes while ensuring safe and efficient training.



REAL-TIME TRACKING OF REPS, SETS, AND PERFORMANCE METRICS

300+ MOVEMENT PROGRAMS

RANGE OF 1 KG TO 75 KG WEIGHT

MULTI-ANGLE ADJUSTABLE BENCH

10" TOUCH SCREEN & WIFI CONNECTED

MULTI-ANGLE ADAPTABILITY

CUSTOMIZABLE ATTACHMENTS

SPACE-EFFICIENT TAKING MINIMAL WALL SPACE

0 / 8  
reps

7.5  
KGS / SIDE

Start

Spotter Calibrate

Eccentric Burnout



- With Aroleap Pro, you can train all major muscle groups through a wide range of strength-focused movements, guided by built-in AI coaching with real-time feedback
- It comes with versatile accessories for bars, grips, and cuffs to support every exercise

# Technical Specifications



## Resistance

2–75 kg adjustable digital resistance



## Accuracy

0.5 kg incremental control



## Motor Technology

BLDC motor (digital weight, no plates)



## Exercise Library

300+ exercises for all major muscle groups



## Smart Training

Eccentric mode, Spotter mode, and AI-guided modes



## Tracking

Reps, resistance, power, and movement quality



## Dimensions

Wall-mounted, compact (~5 ft tall × 2 ft wide)



## Accessories

Barbells, D-handles, ropes, and adjustable arms



## Power Standard

220 V / 5 A home socket



## Digital Control

Adjustment via app or console for precise loading



## Safety

Auto-cutoff to stop resistance instantly during overload



## AI Coaching

Guides workouts, provides feedback, and adapts to progress



We are a complete solution provider for all your Sport Science, Biomechanics, SnC, Physiology, Data Management, Data Analysis, & High Performance Management needs  
We are your partners in excellence!

## Get in Touch...!

**MEDEVICE  
HEALTHTECH**  
TRUST | TECHNOLOGY | INNOVATION

✉ sales@medevice.co

☎ +91 94056 60000 / +91 75881 34567

📍 Pune

## Scan for More Products & Videos



medevice.co



Aroleap Pro Video

f in @ // MedeviceHealthtech