

CABLE

The 1080 Cable is a compact, high-performance training system engineered for strength, speed, and power development. It delivers smooth, controlled resistance with real-time data feedback, enabling athletes and patients to train with precision. Featuring multiple modes—including high-speed, isokinetic, and eccentric overload—it is ideal for both performance training and rehabilitation. Its portable design allows for quick setup on a wall mount, rack, or any training zone, offering a wide range of exercises without occupying excessive space. The result is a versatile platform that elevates the quality and efficiency of every workout.

REAL TIME DISPLAY OF TIME,
DISTANCE, SPEED, FORCE, POWER

7" TOUCH SCREEN &
WIFI CONNECTED

PEDAL CARRYING
HANDLE



DUST & WEATHER
PROTECTED SMART

ANGLE ADJUST
BUTTON $\pm 90^\circ$

3.5 m LINE & LINE
SPEED 5 m/s

GEAR 2 ANCHOR FOR
RESISTANCE UP TO 150 KG

RECHARGEABLE 99
WH BATTERY



- Built-in screen displays current rep, number of reps, best rep, and total work. Measures research grade data in real-time with high-resolution
- The 1080 cable system includes waist belts, harnesses, ankle cuffs, hand grips, bar attachments, pulleys, anchor straps

Technical Specifications



Line Thickness & Length
3.5 mm, 3.5 m
9.8 ft line



Gear 1 Continuous Load
Up to 75 kg (165 lbs)
Max speed 5 m/s (11 mph)



Gear 2 Continuous Load
Up to 150 kg (330 lbs)
Max speed 2.5 m/s (5.5 mph)



Eccentric Overload
Up to 150 kg (330 lbs)



Maximum Speed
5 m/s, 18 km/h, 11 mph



Weight
25 kg (55 lbs)
Without attachments



Dimensions (WxHxD)
370×257×431 mm
14.6×10.1×17.0 in



Data Sampling Rate
200 Hz



Connectivity
Wi-Fi



Touchscreen
7-inch
1300 Lumen



Power
99 Wh Battery
41.4 V, 160 W charger



Port Interface
USB-C



reddot winner 2024



We are a complete solution provider for all your Sport Science, Biomechanics, SnC, Physiology, Data Management, Data Analysis, & High Performance Management needs.
We are your partners in excellence!

Get in Touch...!

**MEDEVICE
HEALTHTECH**

TRUST | TECHNOLOGY | INNOVATION



sales@medevice.co



+91 94056 60000 /+91 75881 34567



Pune

Scan for More Products & Videos



medevice.co



Cable Video



// MedeviceHealthtech