



AROLEAP PRO

The Aroleap Pro is a compact, high-precision device for rehabilitation and functional training. It provides adjustable assistance or resistance to support mobility, strength, and neuromuscular recovery. Therapists can apply controlled loading for safe, progressive rehab and therapeutic exercises. Real-time data tracks movement quality, range of motion, and load for accurate progress monitoring. Its versatile design fits easily in clinics, gyms, or home-care environments. The system helps improve recovery outcomes while ensuring safe and efficient training.



- With Aroleap Pro, you can train all major muscle groups through a wide range of strength-focused movements, guided by built-in AI coaching with real-time feedback
- Accessories: Wireless foot pedal, heavy-duty squat plate, adjustable commercial bench, barbells, D-handles, ropes, and ankle straps

Technical Specifications



Resistance

2–75 kg adjustable digital resistance



Accuracy

0.1 kg incremental control



Motor Technology

BLDC motor (digital weight, no plates)



Exercise Library

300+ exercises for all major muscle groups



Smart Training

Eccentric mode, Spotter mode, and AI-guided modes



Tracking

Reps, resistance, power, and movement quality



Dimensions

Wall-mounted, compact (~5 ft tall × 2 ft wide)



Accessories

Barbells, D-handles, ropes, and adjustable arms



Power Standard

220 V / 5 A home socket



Digital Control

Adjustment via app or console for precise loading



Safety

Auto-cutoff to stop resistance instantly during overload



AI Coaching

Guides workouts, provides feedback, and adapts to progress



We are a complete solution provider for all your Sport Science, Biomechanics, SnC, Physiology, Data Management, Data Analysis, & High Performance Management needs
We are your partners in excellence!

Get in Touch...!

**MEDEVICE
HEALTHTECH**
TRUST | TECHNOLOGY | INNOVATION

✉ sales@medevice.co
☎ +91 94056 60000 / +91 75881 34567
📍 Pune

Scan for More Products & Videos



medevice.co



Aroleap Pro Video

[f](#) [in](#) [@](#) // MedeviceHealthtech