

EASYANGLE

The Meloq EasyAngle is a digital range-of-motion assessment tool that captures joint movement with clinical precision. It uses an advanced inclinometer to instantly measure angles and provide real-time feedback. Therapists and clinicians use it to confidently track mobility, limitations, and progress over time. A core USP is its versatility, allowing assessment of all major joints and motion planes (flexion, extension, rotation, side-bending). Suitable for both upper and lower limbs, EasyAngle delivers repeatable, objective mobility data. This supports individualized treatment planning and informs critical return-to-sport decisions.

ABLE TO MEASURE ALL JOINTS
IN ALL PLANES OF MOTION

SAVES THE LAST 5
MEASUREMENTS

MEASURES FROM 0 TO 180
DEGREES RANGE OF MOTION

DESIGNED FOR PHYSICAL THERAPISTS,
SPORTS CLINICIANS, REHAB SPECIALISTS

REAL-TIME DISPLAY & FEEDBACK ANGLE
SHOWN IMMEDIATELY ON SCREEN

2 X ALIGNMENT GUIDES,
280MM AND 195 MM

DESIGNED FOR ONE
HAND OPERATION

LIGHTWEIGHT, EASY TO CARRY,
SUITABLE FOR CLINIC, FIELD, GYM



- Optimized for assessing multi-joint movements, EasyAngle is trusted by performance coaches, sport therapists, and trainers
- Its supported by straps, anchors, and mounting accessories for accurate testing

Technical Specifications



Operating Environment

Not for home use



Operating Temperature

0°C to 45°C



Enclosure

IP4X



Operating Time

2 weeks



Stand-by Time

11 weeks



Charging Time

2 hours



Sensor Accuracy

±1° within 180°



Certification

Class 1 medical device —
MDD 93/42/EEC (LVFS
2003:11)



We are a complete solution provider for all your Sport Science, Biomechanics, SnC, Physiology, Data Management, Data Analysis, & High Performance Management needs. We are your partners in excellence!

Get in Touch...!

**MEDEVICE
HEALTHTECH**

TRUST | TECHNOLOGY | INNOVATION

✉ sales@medevice.co

☎ +91 94056 60000 / +91 75881 34567

📍 Pune

Scan for More Products & Videos



medevice.co



EasyAngle Video



// MedeviceHealthtech