



picoblu
small • reliable • connected

HEADQUARTERS

Cometa srl
Via Gran Sasso, 18
20010 Bareggio (MI)
Phone: +39 0291410582
Email: sales@cometasystems.com

US BRANCH

Cometa Systems, Inc.
21102 County Road 7560
Newburg, MO 65550-9320

Brochure version 2.0
15 05 2024

cometa

small • reliable • connected

Clip charging

Brand new charger with USB-C connection.
Get a full charge in 3 hours.



18 hours of battery life

8 grams

0 wires



Technical specifications

- The smallest sensor in the world, less than 8 grams
- Bluetooth connectivity
- Entirely without cables
- Compatible with our well tested and reliable EMG & Motion Tools
- Easily variable intra-electrode distance with specific electrodes
- Extremely low artifacts thanks to the symmetrical construction
- More than 18 hours of battery life



Software development kit

Seamless pairing of up to 4 sensors with your laptop, tablet or other device. With software development kit for custom applications.